




## Mixers, Rotators and Rollers

Our rotators give gentle but effective mixing of biological samples to be held in 1.5mL to 50mL micro tubes, used in a variety of applications including immune precipitations, prevention of blood coagulation, latex diagnostics etc.



**MX-RL-E**  
Classic Rotator

**Features:**

- Give gentle but effective mixing of biological samples
- Compact space-saving design fits easily in the incubator
- End-over-end or rolling action, gives gentle but effective mixing
- Used with a wide variety of accessories, for micro tubes with a volume range of 1.5mL-50mL
- Adjustable speed range of 0 to 80 rpm



**MX-RL-Pro**  
LCD Digital Rotator

**Features:**

- The rotator comes with a rotisserie to be holder
- Electronic adjustment of speed and timer
- LCD displays speed and time adjustment
- Speed range from 10 to 70rpm
- End-over-end action, gives gentle but effective mixing
- Used with a wide variety of accessories



	MX-RL-Pro	MX-RL-E
Voltage [VAC]	100-240	100-240
Frequency [Hz]	50/60	50/60
Power [W]	40	20
Motor type	DC motor	DC motor
Tilt angle [o]	0-90	90
Speed range [rpm]	10-70	0-80
Speed display	LCD	-
Timer	Yes	-
Timer display	LCD	-
Time setting range [min]	1-1199	-
Operation type	Continuous / timed operation	Continuous
Dimensions [DxWxH mm]	220x510x260	150x530x190
Weight [kg]	8	1.7
Permissible ambient temperature [°C]	5-40	5-40
Permissible relative humidity	80%	80%
Protection class acc. To DIN EN60529	IP21	IP21

Line ID: @neonics



ติดต่อบริษัท นีโอนิกส์ จำกัด  
Tel: 098-479-5684 หรือ 061-8268939  
E-mail: sale@tools.in.th หรือ sale@neonics.co.th



## Accessories



**Cat. No:18900142**

Used for MX-RL-Pro Rotisserie accessory, for 1.5mL x 48 centrifuge tubes held horizontally



**Cat. No:18900143**

Used for MX-RL-Pro Rotisserie accessory, for 15mL x 24 centrifuge tubes held horizontally



**Cat. No:18900144**

Used for MX-RL-Pro Rotisserie accessory, for 50mL x 24 centrifuge tubes held horizontally



**Cat. No:18900145**

Used for MX-RL-Pro Rotisserie accessory, for 1.5mL x 32 centrifuge tubes held vertically



**Cat. No:18900146**

Used for MX-RL-Pro Rotisserie accessory, for 15mL x 16 centrifuge tubes held vertically



**Cat. No:18900147**

Used for MX-RL-Pro Rotisserie accessory, for 50mL x 16 centrifuge tubes held vertically



Line ID: @neonics



ติดต่อบริษัท เนโอเนิกส์ จำกัด  
Tel: 098-479-5684 หรือ 061-8268939  
E-mail: sale@tools.in.th หรือ sale@neonics.co.th

www.tools.in.th